

Indel Calling Pipeline in the GATK

Guillermo del Angel, Ph.D.

Genome Sequencing and Analysis Group
Medical and Population Genetics Program

Feb 17, 2011



BROAD
INSTITUTE

What are the GATK's indel processing abilities?

GATK Tool	Function
IndelRealigner	Runs multiple sequence alignment on reads and forms consensus indels suitable for variant genotyping.
UnifiedGenotyper	Determines consensus alternate alleles, optimal allele frequency distribution, determines whether sites should be called, assigns genotypes and annotations.
VariantFiltration	Filters calls based on given expressions.
VariantEval	Indel metrics and stratifications for analysis


Step 1: BAM data processing



- Indel realignment is a critical step in preparing BAM's for indel calling.
- We recommend full indel realigning (Smith Waterman) at all sites, **realignment using only known sites is not enough!**

Note: Exome BAM's coming out of Picard have already been fully indel-realigned!

Step 2: Indel discovery



- The genotype likelihoods calculation is inspired by Dindel (with kind permission from C Albers and R Durbin).
- Typical command line:


```
java -jar GenomeAnalysisTK.jar -R ref.fasta -T  
UnifiedGenotyper -L mytargets.list -I myreads.bam  
-o mycalls.vcf -B:dbsnp,VCF dbsnp.vcf -glm DINDEL
```

Only difference with SNP calling!

Some details and caveats...

- All standard parameters used in UG for SNP calling are also valid for indels!
 - E.G. `–stand_call_conf` for a calling threshold.
- Heuristic for controlling sensitivity:
 - We'll only consider indels for genotyping if they are present in N reads, controlled by `–minIndelCnt` parameter. Default value: 5, may want lower value for higher sensitivity in lowpass samples.
- Limitations:
 - Only bi-allelic sites considered. If more than 2 alt. alleles detected at a site, the one with most supporting reads taken.
- NOTE: Application of BAQ will severely degrade indel caller performance. Make sure argument `–baq` is either not included or set to OFF!


Step 3: variant filtration (indels)



- Hard filters are needed for eliminating calls coming from read artifacts.
- This is an ongoing area of improvement, stay tuned on the GATK Wiki for best practice recommendations!
- Example command line with current best practice:

```
java -jar ./dist/GenomeAnalysisTK.jar -T VariantFiltration ref.fasta -o out.vcf -B:variant,VCF input.vcf \
--filterExpression "QUAL<30.0" --filterName "LowQual" \
--filterExpression "SB>=-1.0" --filterName "StrandBias" \
--filterExpression "QD<1.0" --filterName "QualByDepth" \
--filterExpression "(MQ0 >= 4 && ((MQ0 / (1.0 * DP)) > 0.1))" --filterName "HARD_TO_VALIDATE" \
--filterExpression "HRun>=15" --filterName "HomopolymerRun"
```

Step 4 (Optional): Genotype refinement



- Beagle can be used to refine genotypes of indel calls. Current recommended best practice is to merge Indel and SNP calls and running Beagle on combined set. More details our Wiki page.

Assessing indel callsets

- How do we know if the callset that we have is of high sensitivity and high specificity?
- How many variants should we typically get?
- How should indels be distributed in size, allele frequency and types of indels?

VariantEval's support for Indels

```
java -jar GenomeAnalysisTK.jar -B:eval,V
CF mycalls.vcf -T VariantEval -R reffile.fasta -EV
IndelMetricsByAC -EV IndelStatistics -B:dbsnp,VCF
dbsnp.vcf -o output.txt
```

Key module! Produces indel size distributions as well as classification tables

This produces a GATK report file with aggregated statistics.

```
##:GATKReport.v0.1 CountVariants : Counts different classes of variants in the sample
CountVariants CompRod CpG EvalRod JexlExpression Novelty nProcessedLoci nCalledLoci nRefLoci nVariantLoci vari
antRate variantRatePerBp nSNPs nInsertions nDeletions nComplex nNoCalls nHets nHomRef nHomVar nSingleton
s heterozygosity heterozygosityPerBp hetHomRatio indelRate indelRatePerBp deletionInsertionR
atio
CountVariants dbsnp CpG eval none all 63025520 1215 0 1215 0.00
001928 51872.00000000 0 611 604 0 0 724 0 491 0
0.00001149 87051.00000000 1.47454175 0.00001928 51872.00000000 0.98854337

CountVariants dbsnp CpG eval none known 63025520 1000 0 1000 0.00
001587 63025.00000000 0 491 509 0 0 567 0 433 0
0.00000900 111156.00000000 1.30946882 0.00001587 63025.00000000 1.03665988

CountVariants dbsnp CpG eval none novel 63025520 215 0 215 0.00
000341 293141.00000000 0 120 95 0 0 157 0 58 0
0.00000249 401436.00000000 2.70689655 0.00000341 293141.00000000 0.79166667

CountVariants dbsnp all eval none all 63025520 13580 0 13580 0.00
021547 4641.00000000 0 6649 6931 0 0 8852 0 4728 0
0.00014045 7119.00000000 1.87225042 0.00021547 4641.00000000 1.04241239
```


How many indels should I get?

Published estimates In Whole Genome ~ 1 indel/8000 bp

Empirical exome estimate: ~500 indels/exome (33 Mbp)

Lowpass example AC distribution (Chr20)


Total indels by Allele Count, pop= ASN, N=266, 1/Het=8000.0



High heterozygosity in lowpass call set is leading us to focus on improving specificity of our calls.

Exome example AC


Total indels by Allele Count, target captured exomes, N=96, 1/Het=65916.9



High number of singletons: maybe also a lot of false positives?


A typical plot of indel size distribution in whole genome sets

Indel Size Distribution for low-pass 1000G samples, GATK, pop = ASN




Indel size distribution in exomes

Indel Size Distribution for 96 exome capture 1000G samples, GATK



Different indel types come at different rates

Total indels by Indel type, pop= ASN, N=266



These counts for the Phase 1 1000 Genome samples taken for illustration only

Consistent difference in the number of A,T vs. C,G repeat expansions: Ratio of (A+T)/(C+G) expansions can be used as a specificity metric!

Other callers

- Aside from the GATK, SAMTools and DINDEL can be alternatively used for indel calling.
- Example command line using SAMTools' mpileup caller:

```
samtools mpileup -ugf ref.fasta reads.bam | ../samtools/  
bcftools/bcftools view -vc - > myout.vcf
```

- More info at:

- <http://samtools.sourceforge.net/mpileup.shtml>
- <http://www.sanger.ac.uk/resources/software/dindel/>