



detecting, annotation and visualization of the alternative splicing from RNA-Seq data

Home	Download	Screenshots	Document	Useful links	Feedback	About us
----------------------	--------------------------	-----------------------------	--------------------------	------------------------------	--------------------------	--------------------------

Requirements

- SplicingViewer is a software integrating a command line program, which aims to detect splice junctions and annotate alternative splicing event patterns, and viewer for visualizing the alternative splicing event patterns. The program and viewer are both written in java programming language and require java [version 1.6.0 or higher](#).

Pipeline to detect and annotate the alternative splicing

- [Step 1: Mapping short reads to the reference genome sequence](#)
Firstly, the RNA-Seq reads are mapped to the reference genome sequence with mapping tool by the users, which can generate SAM/BAM format results or whose results can be converted into SAM/BAM format. The user can choose the mapping tool by themselves. The tolerances should be set to allow at most two mismatches. According to the validation between different mapping tools, we recommend users to use [BWA](#) as the aligner, not only because it relative high accuracy and multiple threads abilities, but it can also support native SAM/BAM result generating. The conversion of SAM to BAM, sort and index of BAM file can be achieved by [SAMtools](#). Example shell script using BWA as mapping tool can be found [here](#).
- [Step 2: Getting unmapped short reads](#)
After the mapping of short reads to the reference genome, the mapped records contained in the result BAM file in the previous step are used to calculate the depth of coverage of each locus contained in the RefGene model bound with the GATK (The Genome Analysis Toolkit, incorporated into the AlternativeSplicing.jar) in next step. The unmapped records should be converted into fastq read sequences by [SamToFastq.jar](#) from [Picard](#). Example shell script can be found [here](#).
- [Step 3: Detecting splice junctions](#)
The step 'sp' contained in [AlternativeSplicing.jar](#) were used to get the candidate splice junctions. Input of this step contains reference genome sequence, RefGene in GFF format ([Example](#), generated from human UCSC RefGene with our python [script](#)) or BED format ([Example](#)), the sorted and indexed BAM result generated in step 1. This step will generate three files in the current "SPresult" directory, which are 'spliceSites.txt ([Example](#))', 'spliceJunctions.fa ([Example](#))', 'transShortRun.txt ([Example](#))'. Example script to run this step can be found [here](#).
- [Step 4: Mapping the unmapped short reads to splicing junction sequence](#)
Map the unmapped reads, generated in the step 2, to splice junction fasta sequences (splicingJunctions.fa). The tolerances should be set to allow at most two mismatches. The SAM result of this step will be used for the following step, and need not to be converted to BAM. Example shell script using BWA as mapping tool can be found [here](#).
- [Step 5: Annotating alternative splicing event patterns](#)
The step 'as' contained in [AlternativeSplicing.jar](#) were used to annotate the seven alternative splicing patterns, which are skipped exons (SE), mutually exclusive exons (MXE), alternative 5' splicing sites (A5SS), alternative 3' splicing sites (A3SS), alternative first exons (AFE), alternative last exons (ALE), and retained intron (RI). The input of this step contain four files: the RefGene file (GFF or BED format, same as the RefGene file used in step 3), the splicing site text file (spliceSites.txt) generated by step 3, the transcription short run file (transShortRun.txt) generated in step 3, and the SAM file generated in step 4 and contain the reads mapping to splicing site junction sequence. This step will generate the AS patterns result file, whose name can be assigned by the user or named 'AS_result.txt' by default, and the indexed bam file which contain the supported junction reads information. User can also get the splice junction bed format file by add -j option in AlternativeSplicing.jar command line. All the results will be contained in the current AS Result directory. Example shell script can be found [here](#).
- [Step 6: Use viewer to visualize the alternative splicing patterns](#)
The input of the viewer contains five files: the reference genome sequence, the sorted and indexed BAM file generated in step 1, the sorted and indexed BAM file containing the junction reads supporting the AS result, the RefGene in GFF or BED format, and the AS result file generated by previous step.

Overview introduction of the viewer

- [Menus](#)

File	
New Session	
Load Session	Ctrl-O
Save Session	Ctrl-S
Save Session as	Ctrl+Shift-S
Export AS Data	
Export Annotation Image	
Export Alignment Image	
Close	
Quit	Ctrl-Q

New session – for first time use, click new session to provide the following files:

1. Reference - the reference genome sequence in fasta format. Genome index file <ref.fa>.fai and directory file <ref>.dict should be supplied in the same directory of reference genome, which can be created by 'samtools faidx' and [CreateSequenceDictionary.jar](#) from [Picard](#), respectively. Dictionary file can also be created by the viewer automatically.
2. RG Alignment - a alignment file in BAM format, containing the short reads mapped to the reference genome. The BAM file should be firstly sorted using 'samtools sort', then be indexed using 'samtools index', resulting index file <aln.bam>.bai in the same directory of BAM file.
3. SJ Alignment - an alignment file containing the unmapped reads mapping to the splice junction sequences in BAM format. The according SAM result is one of results of command line program, then should be converted, sorted and indexed use [SAMtools](#) , resulting <junctionRecord.bam>.bai in same directory.
4. RefGene - known gene model annotation files in GFF ([example](#)) or BED ([Example](#)) format. Please refer to the example gene model file for the format.
5. ASResult - alternative splicing event patterns annotation result file in GFF format generated by the pipeline.

Click "Display" button, loading file process will be finished in a few seconds. Then the main page will be displayed.

Setting Path

Reference	E:\test\hg18.fa	+
RG Alignment	E:\test\total.sort.bam	+
SJ Alignment	E:\test\JunctionRecord.sort.bam	+
RefGene	E:\test\refGene_exon.gff	+
ASResult	E:\test\AS_events_AI.gff	+

Display

Load session - quickly open an existing session saved as a log file.

Save session - save current session.

Save session as - allows you to choose a new name, directory to save the session

Export AS data- export the AS pattern result as the tabulated format.

Export Annotation image- export annotation track image in PDF format.

Export Alignment image- export alignment track image in PDF format.

Close - close current session.

Quit - exit the application.

• [View Menu](#)

View

- ☒ Show Annotation
- ☒ Show Coverage
- ☒ Show Reference
- ☒ Show Alignment

Showing annotation: option for showing annotation track and gene and AS lists.

Showing Coverage: option for showing coverage

Showing Reference: option for showing reference

Showing Alignment: option for showing alignment

• [Settings Menu](#)

Settings

Annotation Parameters

Alignment Parameters

Project Parameters

Annotation parameters: provide user-defined parameters highly configurable allowing users to define view mode, arrow mode, vertical space, change track display height and background color. Annotation tracks: functions for reordering of track location and adding/deleting related tracks.

Annotation Option

View

View Mode

Nooverlap

Single Line

Arrow Mode

Show Gene Arrow

Hide Gene Arrow

Track view

Show GENE

Show AS

Track show order

Gene Track up

As Track up

Unit Height

Track Height

Unit Vertical Space

Track Background

Alignment parameters: users can change font and colors in many different combinations, such as each base and background color. Zooming can be performed by choosing from the settings menu or using toolbar buttons "+" and "-".

Alignment Option

Font

FamilyDialog

Size

Color

Style

Bold

Italic

Show Base String

Show Reference Base String

Show Simple Read Base String

Show Highlight Read Base String

Show Arrow

Show Border

Read Base Color

Reference And HightLight Base Color

Project Parameters: users can change the quality type, adjust and view windows size, and can set whether to show the pop tip.

Main Option

Global

Quality Type

Illumina

Adjust Window Size


100

View Window Size

1,200

Show PopTip


- [Help](#)




- **Toolbars**
Several basic items can be found on the toolbars. You can use these tool bars for zooming in/out, selecting current chromosome or coordinate settings etc.



- **RefGene id and AS event id panel**
When you have established or opened a session, the RefGene id and AS event id will be shown in the left panel of viewer. The id item can be clicked and the alignment track and annotation track will be jump to the according sites in the right panel. You can also search the gene by the id or name.



- **Annotation track panel**
The AS events and RefGene are shown graphically in this annotation track panel located in up-right of the viewer. you can view the junctions mapping reads when you mouse suspend on the polyline, indicating the junction, then click it when the color of polyline turn to red, the support junction mapping reads will be shown in the popup dialog. The seven alternative splicing event patterns are shown in the following shapes.



- [Alignment track panel](#)
The alignment track panel contains three sub panels: the depth of coverage panel, the reference genome sequence panel, and the alignment panel.
A step by step documentation for the test data application will come soon!